

# Oligonucleotide Analysis by LC-MS/MS

CASE STUDY



## The Oligonucleotides Analysed RESO



Oligonucleotide 1 (Analyte): MW: 6258 Da 18-mer, fully phosphorothioated, 2'-O-methylated modified bases Oligonucleotide 2 (Internal Standard): MW:6015.8 Da 19-mer, fully phosphorothioated

Backbone modification by the incorporation of phosphorothioated linkages was selected as this a very common oligonucleotide modification; it confers nuclease resistance and enhances in vivo stability. The modification results in the addition of a chiral centre presenting chromatographic challenges. Extraction challenges due to being highly protein bound.

Second generation antisense oligonucleotides have a ribose sugar modification to improve plasma stability and increase tissue half-lives. Therefore a 2'-O-methyl modification was elected to reflect the generational evolution in synthetic oligonucleotide therapeutics.

#### **Overcoming Challenges of Oligonucleotide Analysis**





•Polypropylene containers, including mobile phase bottles used throughout. No glassware used. Minimising loss due to absorption and reducing Na adducts in the mass spectrometer.

•Serial dilution in plasma from stock solution performed for preparation of calibration standards and QC samples. Diluted precipitated matrix used to prepare post-spike solutions to overcome non-specific binding.

Extraction Recovery

Optimisation of all steps of the SPE extraction resulting in 70% recovery for both oligonucleotides allowing for a low LLOQ to be achieved.
Addition of the cysteine to the Lysis-Loading buffer to prevent desulfurization during sample processing.

Passivation of
 Initial poor pe
 solvent.
 Mabila phases

Passivation of the LC system and analytical column with ion-pair mobile phase to ensure rapid system equilibration.
Initial poor peak shape due to the phosphorothioate stereocenter overcome by inclusion of 400 mM HFIP in both LC mobile phases and weak wash solvent.

•Mobile phases prepared daily to maintain peak shape and intensity.

Selectivity

•Selectivity can be challenging due to production of small, non-selective product ions after fragmentation by CID. Numerous charge states and product ions were monitored throughout the development allowing for selection of the optimal MRM to achieve selectivity and sensitivity.



#### Methodology



500 µL human plasma



Phenomenex Clarity lysis loading buffer with cysteine



Waters Oasis WAX SPE using

Otto Positive Pressure Manifold





Waters Acquity UHPLC coupled to a Sciex 5500 Triple-Quad MS



### **Performance Data**



Plasma stability demonstrated after:

- 17 hours at room temperature
- 13 days at -80°C
- 1 freeze-thaw cycle

Additional assessments:

- Selectivity in individuals free of interference
- Matrix effects in individuals acceptable



Replicate	QC 1ª	QC 2ª	QC 3 <sup>a</sup>	QC 4 <sup>b</sup>	QC 5 <sup>⊳</sup>	<b>QC 6</b> <sup>b</sup>	<b>Q</b> C 7 <sup>ь</sup>
Mean	0.3	0.5	1.0	4.76	15.6	787	1770
% <b>CV</b>	5.66	3.76	2.68	3.3	4.9	2.0	2.0
Accuracy (%)	0.67	2.4	0.2	-4.83	3.74	-1.59	-5.57

Table 1: Quality Control Precision and Accuracy

<sup>a</sup> Analytical assay 500 μL sample volume, 0.25-250 ng/mL <sup>b</sup> Analytical assay 200 μL sample volume, 5-2500 ng/mL



**Contact us** to discuss how we can support your cutting-edge oligonucleotide therapeutics.

www.resolian.com